

HYDRION® PLASTIC pH STRIPS

Hydrion brand Plastic pH Strips provide instant, easy-to-read results, with clear, bright single color matches across the entire pH range.

Hydrion Plastic pH Strips are packaged in convenient flip-top vials, with matching color charts. Strips are available in 5 pH ranges. Short-range strips offer clear, bright, single-color matches at every 0.2/0.3 interval (5.5–8.0) or every 0.5 interval (0.0–6.0, 5.0–9.0, 6.5–13.0). The wide-range strip provides accurate pH measurement at each 1.0 interval over the entire pH range (0.0–14.0).

- **Fast, accurate, easy-to-use: one dip—one color match**
- **Packaged in convenient flip-top vials**
- **Made in the U.S.A.**



HYDRION® PLASTIC pH STRIPS

- Strip dimensions: 2 3/4 in. x 1/4 in.
- Packaged in flip-top vials
- 100 strips per vial

Catalog No.	pH Range	pH Readings
9200	0.0–6.0	0.0 0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0 5.5 6.0
9400	5.0–9.0	5.0 5.5 6.0 6.5 7.0 7.5 8.0 8.5 9.0
9700	5.5–8.0	5.5 5.8 6.0 6.2 6.4 6.6 7.0 7.2 7.4 7.6 8.0
9600	6.5–13.0	6.5 7.0 7.5 8.0 8.5 9.0 9.5 10.0 10.5 11.0 11.5 12.0 12.5 13.0
9800	0.0–14.0	0.0 1.0 2.0 3.0 4.0 5.0 6.0 7.0 8.0 9.0 10.0 11.0 12.0 13.0 14.0

Sale Unit	Price	Minimum Order Quantity
Case of 6 vials	Distributor pricing available	1 case

How to Measure pH with Hydrion Strips

- Dip the strip into the test solution
- Remove the strip and immediately compare the resulting color with the matching pH color chart



© 2013 Micro Essential Laboratory Inc. All rights reserved. Printed in U.S.A.

For more information contact
 Hawaii Chemical & Scientific at 808-841-4265
Sales@HawaiiScientific.com
www.HawaiiScientific.com